





DIRECT DETECTION OF DNA METHYLATION AT CpA UNITS



Marie W. Johannsen^a, Simon R. Gerrard^b, Tracy Melvin^c, Tom Brown^a, ^aDepartment of Chemistry, University of Oxford, Oxford OX1 3TA, marie.johannsen@chem.ox.ac.uk, ^bSchool of Chemistry and ^cOptoelectronics Research Centre, University of Southampton, Highfield, Southampton SO17 1BJ, UK,

DNA **methylation patterns** are emerging as targets of great biological relevance. We have been working to develop a method to detect methylation patterns directly, without the need for bisulphite conversion or polymerase-aided replication.

PhP, a monomer first investigated for its ability to stabilise triplexes with cytosines in the purine strand (**CG inversions**) has been found to decrease stability of triplexes if the purine strand C is methylated.




While **CpA methylation** is easily detected by this method, it is not possible to distinguish between CpG and ^{Me}CpG. While the methylation frequency of CpG sites is usually the highest,

CpA methylation is common in **embryonic stem cells** and in **plant cells**. In rice, the methylation frequency of CHG sites is 21%, while in maize it is as high as 74%

In conclusion, we have developed a method capable of detecting **CpA methylation** in a purine stretch of a duplex. By using triplex forming oligonucleotides containing the monomer **PhP**, methylation at a CpA cytosine leads to a detectable change in **denaturation temperature**

Triplex DNA

Triplex DNA is formed of **three strands of DNA**. Two strands form a normal duplex with Watson-Crick base pairing, with a third strand bound in the major groove *via* **Hoogsteen** base pairing. It is easiest to form a triplex where the third strand is a polypyrimidine binding to a stretch of purines in the DNA sequence. Binding a triplex-forming oligonucleotide does not require denaturation of the duplex and the **denaturation temperature** (T_m) of the triplex can be measured separately from the duplex denaturation.



References:

- M. W. Johannsen, S. R. Gerrard, T. Melvin, T. Brown, *Chem. Commun.*, 2014, **50**, 551-553
 The ENCODE Project Consortium, *Nature*, 2012, 489, 57-74
 J. Shimooka *et al.*, *Biochemistry*, 2013, **52**, 1344-1353
 M. Frommer *et al.*, *Proc. Natl. Acad. Sci. USA*, 1992, **89**, 1827-1831
 J. I. Gent *et al.*, *Genome Res.*, 2013, **23**, 628-637
 J. L. Asensio *et al.*, *Nucleic Acids Res.*, 1998, **26**, 3677-3686
 S. R. Gerrard *et al.*, *Org. Biomol. Chem.*, 2010, **8**, 5087-5096

5' - CTM CTX TMT MCC TMT MTC
 5' - GGA GAG GAY AGA GGG AGA GAG GGC
 3' - CCT CTC CTG TCT CCC TCT CTC CCG

Example oligonucleotide set.
 M = MeC, X = PhP, Y = C or MeC